

A GloMax[®] 20/20 Luminometer Method for Dual-Luciferase[®] Reporter Assay

INTRODUCTION

The GloMax[®] 20/20 Luminometer in combination with the Dual-Luciferase[®] Reporter (DLR[™]) Assay System provides a convenient, rapid and sensitive procedure for quantifying gene expression. Transcriptional regulation, coupled to the expression of a luciferase reporter gene, is regularly used to study a wide range of biological events in cultured cells. Luciferase is an ideal reporter because of the absence of endogenous luciferase activity in mammalian cells, and the functional enzyme is created immediately upon translation (1,2).

The DLR[™] Assay System contains two different luciferase reporter enzymes that are expressed simultaneously in each cell. Typically, the experimental reporter is correlated with the effect of specific experimental conditions, while the activity of the co-transfected "control" reporter gene provides an internal control, which serves as the baseline response. Normalizing the experimental reporter gene to the activity of an internal control minimizes the variability caused by differences in cell viability and transfection efficiency. Thus, dual reporter assays allow more reliable interpretation of the experimental data by reducing extraneous influences.

The experimental and control luciferase enzymes used in the DLR[™] Assay System have distinct evolutionary origins. The firefly luciferase and the *Renilla* (sea pansy) luciferase can discriminate between their respective bioluminescent substrates and do not cross-activate.

The firefly and *Renilla* substrates have been developed specifically to maximize the sensitivity of the assay reagent. This system is widely used in life science research because of the superior light generation and high signal to noise ratio. Dual-Luciferase[®] Reporter reagents are compatible with commonly used culture media for mammalian cells including RPMI 1640, MEM α , DMEM and Ham's F12.

These reagents are designed for use with the Passive Lysis Buffer that comes with the kit and also is available separately Cat.# E1941.

The sensitivity and wide dynamic range offered by the GloMax[®] 20/20 Luminometer make it highly suited for the DLR[™] Assay application. The GloMax[®] 20/20 software facilitates a quick setup for the DLR[™] Assay with a pre-installed template. The internal automatic injectors are easy to use and completely accessible.

The GloMax[®] 20/20 Luminometer can detect as little as 1×10^{-21} moles of firefly luciferase enzyme using Luciferase Assay Reagent II (LAR II). All tests were conducted using purified recombinant firefly luciferase enzyme (Cat.# E1701) and purified *Renilla* recombinant enzyme (Chemicon Cat.# 4400).




Figure 1. DLR[™] Assay was performed on the GloMax[®] 20/20 Luminometer using the DLR Assay System with recombinant firefly luciferase (1×10^{-21} moles to 1×10^{-14} moles) and *Renilla* luciferase (1×10^{-16} moles).

MATERIALS REQUIRED

- GloMax[®] 20/20 Luminometer
- Microcentrifuge Tube Holder
- 1.5 mL Microcentrifuge Tubes
- Dual-Luciferase[®] Reporter Assay System (Cat.# E1980)
- p200 pipette and pipette tips
- p20 pipette and pipette tips

EXPERIMENT PROTOCOL

1. Reagent Preparation

Luciferase Assay Buffer II and Luciferase Assay Substrate:

Use as supplied. Store at –20°C, where it is stable for up to 6 months. The Luciferase Assay Substrate also may be stored at 4°C for up to 1 month.

Transfer the contents of one bottle of Luciferase Assay Buffer II into one vial of Luciferase Assay Substrate. Mix by inversion until the substrate is thoroughly dissolved. Use reconstituted Luciferase Assay II Reagent (LAR II) on the same day it is prepared or divide into working aliquots and store at –20°C for 1 month or at –70°C for up to 1 year.

Stop & Glo[®] Substrate and Stop & Glo[®] Buffer: Use as supplied. Store below –20°C.

Stop & Glo[®] Buffer Substrate Solvent: Use as supplied. Store below 25°C.

To make **Stop & Glo[®] Reagent**, dilute the 50X Stop & Glo[®] Substrate to 1X concentration using Stop & Glo[®] Buffer in a glass or siliconized polypropylene tube. Mix by inversion. Use reconstituted Stop & Glo[®] Substrate on the same day it is prepared or store at –20°C for up to 2 weeks.

Passive Lysis Buffer: To make 1X Passive Lysis Buffer, dilute the 5X Passive Lysis Buffer in deionized water. Store below 25°C.

Note: The temperature of the Luciferase Assay Buffer II and Stop & Glo[®] Buffer should be held constant at room temperature while quantifying luminescence, since luciferase activity is temperature-dependent. Reagent stored frozen after reconstitution must be thawed below 25°C to ensure reagent performance. Mix well after thawing. The simplest method for thawing is placing the reagent in a water bath at room temperature.

2. Instrument Setup

- Turn ON the GloMax[®] 20/20. A 5-minute warm-up period is recommended but not necessary.
- Touch "Run Promega Protocol" from the "Protocols" menu.
- Select "DLR-2-INJ" from the list of Promega protocols. The "Parameters" screen appears next with pre-programmed settings that are optimized for the DLR[™] Assay with two automatic injectors.
- Touch "OK" to go to the "Home" screen.

3. Sample Analysis

- Remove the cell cultures from the incubator. **Note:** For maximum reproducibility, equilibrate cell cultures to room temperature before adding reagent.
- Prepare the injectors. Place the inlet tubing for injector 1 into the bottle of LAR II. Place the inlet tubing for injector 2 into the bottle of Stop & Glo[®] Reagent. Place a waste container underneath the injector tips.
- Touch "Inj Func" and prime both injectors. **Note:** Do not switch injectors. Residual Stop & Glo[®] Reagent will quench the firefly luciferase reporter activity. We recommend dedicating injector 1 for LAR II and injector 2 for Stop & Glo[®] Reagent.
- Add 20 µL of your sample into a microcentrifuge tube.
- Place the microcentrifuge tube into the tube holder, and close the lid.
- Touch "Measure Luminescence" to begin measurement. The GloMax[®] 20/20 will automatically inject 100 µL of LAR II into the sample. After a 2-second delay, the GloMax[®] 20/20 will measure the firefly signal over a 10-second integration period. Then, the GloMax[®] 20/20 will automatically inject 100 µL of Stop & Glo[®] Reagent into the microcentrifuge tube. After a 2-second delay, the GloMax[®] 20/20 will measure the *Renilla* signal over a 10-second integration period.
- Repeat previous step for each sample.
- When all measurements are complete, flush the injectors with deionized water followed by 70% ethanol. Perform a final flush with deionized water and allow the water to remain in the injector system.

REFERENCES

1. Ow, D.W. *et al.* (1986) Transient and stable expression of the firefly luciferase gene in plant cells and transgenic plants. *Science* **234**, 856–9.
2. De Wet, J.R. *et al.* (1987) Firefly luciferase gene: structure and expression in mammalian cells. *Mol. Cell. Biol.* **7**, 725–37.

DLR, Dual-Luciferase, GloMax and Stop & Glo are trademarks of Promega Corporation.

CAUTION: The lyophilized Luciferase Assay Substrate contains dithiothreitol (DTT) and is therefore classified as hazardous. The reconstituted reagent is not known to present any hazards as the concentration of DTT is less than 1%. However, we recommend the use of gloves, lab coats and eye protection when working with these or any chemical reagents. Promega assumes no liability for damage resulting from handling or contact with these products.

CONTACT INFORMATION

Toll-Free: (800) 356-9526
Fax: (800) 356-1970

www.promega.com

Email: custserv@promega.com

Mailing Address:

Promega Corporation
2800 Woods Hollow Rd.
Madison, WI 53711 USA